

eliminate matrix interferences, the analysis is fairly rapid and reliable, at least 4-5 times more rapid than the accepted atomic absorption method of the AOAC (1973).

#### LITERATURE CITED

- Amos, M. D., *Amer. Lab.*, 4(8), 57 (1972).  
 Analytical Methods Committee, *Analyst (London)* 79, 397 (1954).  
 Associated Press release, *The Cincinnati Enquirer* (Sept 18, 1973).  
 Association of Official Analytical Chemists, *J. Ass. Offic. Anal. Chem.* 56, 481 (1973).  
 Brodie, K. G., Matousek, J. P., *Anal. Chem.* 43, 1557 (1971).  
 Christian, G. D., *Anal. Chem.* 41, 24A (1969).  
*Consumer Reports*, 601 (Oct, 1973).  
 Donega, H. M., Burgess, T. E., *Anal. Chem.* 42, 1521 (1970).  
 Fiorino, J. A., Moffitt, R. A., Woodson, A. L., Gajan, R. J., Huskey, G. E., Scholz, R. G., *J. Ass. Offic. Anal. Chem.* 56, 1246 (1973).  
 Gajan, R. J., Larry, D., *J. Ass. Offic. Anal. Chem.* 55, 727 (1972).  
 Hoover, W. L., Melton, J. R., Howard, P. A., *J. Ass. Offic. Anal. Chem.* 54, 860 (1971).  
 Hwang, J. Y., Mokeler, C. J., Ullucci, P. A., *Anal. Chem.* 44, 2018 (1972).  
 Hwang, J. Y., Ullucci, P. A., Smith, S. B., Jr., Malenfant, A. L., *Anal. Chem.* 43, 1319 (1971).  
 Kurz, D., Roach, J., Eyring, E. J., *Anal. Biochem.* 53, 586 (1973).  
 Manning, D. C., *Amer. Lab.*, 5(8), 37 (1973).  
 Matousek, J. P., *Amer. Lab.*, 3(6), 45 (1971).  
 Matousek, J. P., Stevens, B. J., *Clin. Chem.* 17, 363 (1971).  
 Munns, R. K., Holland, D. C., *J. Ass. Offic. Anal. Chem.* 54, 202 (1971).  
 Shea, K. P., *Environment* 15, 6 (March, 1973).  
 Slavin, W., "Atomic Absorption Spectroscopy," Interscience, New York, N. Y., 1968, p 118.  
 Takeuchi, T., Yanagisawa, M., Suzuki, M., *Talanta* 19, 465 (1972).  
 Taphorn, J. E., III, Department of Chemistry, University of Cincinnati, 1973, private communication.  
 United Press International release, *The Cincinnati Enquirer* (Sept 19, 1973).  
 Yanagisawa, M., Takeuchi, T., Suzuki, M., *Anal. Chim. Acta* 64, 381 (1973).  
 Youden, W. J., "Statistical Methods for Chemists," Wiley, New York, N. Y., 1955, pp 43-45.

Received for review October 31, 1973. Accepted May 22, 1974. The authors gratefully acknowledge the National Science Foundation who supported this study in the form of a research grant (GP36478).

## Volatile Constituents of Pressure Cooked Pork Liver

Cynthia J. Mussinan\* and John P. Walradt

The volatile constituents of pressure cooked pork liver were isolated by simultaneous steam distillation and continuous solvent extraction. Analysis by gas chromatography and coupled gc-mass

spectrometry led to the identification of 179 compounds. The mass spectral identifications were confirmed by matching retention indices. Pyrazines were the largest group of compounds found.

As a food, liver is an important source of protein, fat, and vitamins, particularly vitamins A, D, E, and vitamin B complex. A search of the literature, however, has revealed that no previous work has been reported on the volatile constituents of cooked liver. Pork liver was chosen for this study because its odor is typical of liver, and it is considerably stronger than calf or beef liver.

#### EXPERIMENTAL SECTION

Locally procured pork liver (42.5 lb) was sliced into small pieces and passed through a Fitz Mill Model D Comminuter (knives forward, no screen). The sliced liver was transferred to a 20-gal. steam jacketed, doubly stirred, stainless steel, pressure reaction vessel (Groen Div., Dover Corporation, Elk Grove, Ill.) and then slurried with 27 l. of distilled water. The vessel was sealed, heated to 325°F, and held at this temperature for 15 min. The maximum head pressure attained was 98 psi. At the end of the cooking period the vessel was cooled by passing water through the jacket, and the contents, which had a typical cooked liver aroma, were filtered through cheese cloth and a wire basket funnel into 5-gal polyethylene containers and stored at -20° until used.

The filtrate was atmospherically steam distilled, and the distillate continuously extracted with distilled diethyl ether in a scaled-up (22 l.) model of William's apparatus (Williams, 1969). The distillation-extraction was carried out over a 48-hr period for each of the two batches.

The extracts were combined, dried over anhydrous sodium sulfate, and initially concentrated to about 80 ml by careful distillation in a 1-l. Kuderna-Danish concentrator

(Kontes Glass Co., Vineland, N. J.) equipped with a 508 mm × 25.4 mm i.d. reflux column packed with 6 mm × 6 mm Raschig rings. The extract was further concentrated to about 15 ml in a 100-ml round-bottomed flask with a 35 mm × 10 mm o.d. test tube sealed to the bottom and equipped with a 300 mm × 13 mm i.d. Vigreux reflux column. The open container was then allowed to stand at room temperature until the sample had concentrated to a final volume of about 10 ml. The concentrated extract possessed a typical cooked liver aroma.

The concentrate was analyzed on a Hitachi Model RMU-6E mass spectrometer coupled with a Hewlett-Packard Model 5750 gas chromatograph using a Watson-Biemann helium separator (Watson and Biemann, 1965). The chromatographic columns used were 1000 ft × 0.03 in. i.d. stainless steel open tubular columns coated with SF-96 and Carbowax 20M. Further analysis was preceded by area trapping in Varian 1-ml collection bottles cooled with Dry Ice-isopropyl alcohol. The instrument used was a Varian Series 712 preparative gas chromatograph with a flame ionization detector employing a 12 ft × 3/8 in. stainless steel column packed with 20% SE-52 on 45-60 mesh acid-washed, DMCS-treated Chromosorb W. The oven temperature was programmed from 80 to 225° at 4°/min after a 5-min post-injection hold. The injector and detector temperatures were 230° and the helium flow rate was approximately 300 ml/min. Six traps were collected at arbitrary intervals and analyzed in the gc-mass spectrometry system. Traps 1-4 were analyzed on 1000 ft × 0.03 in. stainless steel open tubular SF-96 and Carbowax 20M columns. Trap 5 was analyzed on a 500 ft × 0.03 in. stainless steel open tubular SF-96 column and an 8 ft × 1/8 in. stainless steel column packed with 10% Carbowax 20M on 80-100 mesh acid-washed, DMCS-treated Chromosorb W.

\*International Flavors and Fragrances, Inc., Union Beach, New Jersey 07735.

Table I. Compounds Identified from Cooked Pork Liver

| Compound                    | Retention index <sup>a</sup> |         |       |         | Mass spectral data <sup>b</sup> |
|-----------------------------|------------------------------|---------|-------|---------|---------------------------------|
|                             | CBW                          |         | SF-96 |         |                                 |
|                             | Known                        | Unknown | Known | Unknown |                                 |
| Hydrocarbons                |                              |         |       |         |                                 |
| 2-Pentene                   |                              | <1.00   |       |         | 55, 70, 42, 41, 39, 29          |
| Heptane                     |                              |         | 3.05  | 3.05    | 43, 41, 29, 27, 57, 100         |
| Pentadecene                 |                              |         |       | 10.77   | 41, 55, 57, 69, 83, 56... 210   |
| Pentadecane                 |                              | 9.10    |       |         | 57, 71, 43, 41, 85, 55... 212   |
| Hexadecene                  |                              | 10.43   |       |         | 41, 43, 55, 57, 83, 69... 224   |
| Hexadecane                  |                              |         |       |         | 57, 43, 71, 41, 29, 85... 226   |
| Heptadecane                 |                              |         |       | 12.38   | 57, 43, 71, 41, 85, 55... 240   |
| Benzene                     | 2.92                         | 2.80    |       |         | 78, 51, 52, 50, 77              |
| Ethylbenzene                |                              |         | 4.74  | 4.76    | 91, 106, 51, 39, 65, 27         |
| Hexylbenzene                |                              |         |       |         | 43, 77, 91, 162, 27, 29         |
| Toluene                     |                              |         | 3.75  | 3.62    | 91, 92, 39, 65, 61              |
| <i>m</i> -Xylene            | 5.17                         | 5.17    |       |         | 91, 106, 105, 39, 51, 77        |
| Limonene                    |                              |         | 6.49  | 6.33    | 68, 93, 39, 67, 41, 27... 136   |
| Alcohols                    |                              |         |       |         |                                 |
| Methanol                    |                              |         | <1.00 | <1.00   | 31, 32, 29, 28, 18              |
| Ethanol                     |                              |         | <1.00 | <1.00   | 31, 45, 29, 27, 46              |
| 2-Propanol                  | 2.30                         | 2.26    |       |         | 45, 43, 27, 29, 41... 60        |
| 1-Butanol                   | 5.00                         | 5.00    |       |         | 56, 31, 41, 43, 27, 42... 74    |
| 2-Methyl-1-propanol         | 4.37                         | 4.38    |       |         | 43, 31, 41, 42, 27, 33... 74    |
| 2-Butanol                   | 3.55                         | 3.66    |       |         | 45, 27, 59, 31, 29, 43... 74    |
| 3-Methyl-1-butanol          | 5.55                         | 5.50    |       |         | 41, 29, 55, 31, 57... 88        |
| 2-Methyl-2-butanol          | 3.45                         | 3.40    |       |         | 59, 73, 43, 55, 31... 88        |
| 2-Methyl-3-buten-2-ol       | 3.82                         | 3.85    |       |         | 71, 43, 59, 41, 39... 86        |
| 1-Octen-3-ol                | 7.97                         | 8.00    |       |         | 57, 43, 41, 72, 55, 27... 128   |
| Aldehydes                   |                              |         |       |         |                                 |
| Acetaldehyde                |                              |         | <1.00 | <1.00   | 29, 44, 43, 26, 42, 27          |
| Propanal                    | <1.00                        | <1.00   |       |         | 29, 28, 27, 58, 26, 57          |
| Butanal                     |                              |         | 1.79  | 1.76    | 44, 43, 72, 41, 27, 29          |
| Pentanal                    |                              |         | 2.79  | 2.81    | 44, 29, 27, 41, 58, 28... 86    |
| Hexanal                     |                              |         | 4.00  | 4.00    | 41, 29, 43, 27, 56, 57... 100   |
| Heptanal                    |                              |         | 5.11  | 5.00    | 29, 27, 41, 44, 43... 114       |
| Octanal                     | 6.59                         | 6.67    |       |         | 43, 29, 41, 44, 57, 55... 128   |
| Nonanal                     |                              |         | 7.20  | 7.17    | 41, 43, 57, 56, 44, 29... 142   |
| Dodecanal                   |                              |         | 10.14 | 10.26   | 43, 41, 57, 82, 29, 68... 184   |
| Tridecanal                  | 11.78                        | 11.70   |       |         | 57, 43, 41, 82, 55... 198       |
| Tetradecanal                |                              |         | 12.26 | 12.13   | 43, 41, 57, 55, 29, 82... 212   |
| Pentadecanal                |                              |         |       |         | 43, 41, 57, 55, 29, 82... 226   |
| Hexadecanal                 |                              |         |       |         | 43, 41, 57, 55, 82, 29... 240   |
| Heptadecanal                |                              |         |       |         | 43, 41, 57, 55, 29, 82... 254   |
| Octadecanal                 |                              |         |       |         | 43, 41, 57, 55, 82, 29... 268   |
| Isobutyraldehyde            |                              |         | 1.28  | 1.31    | 43, 41, 72, 27, 29              |
| Isovaleraldehyde            |                              |         | 2.37  | 2.35    | 44, 41, 43, 46, 58, 29... 86    |
| 2-Methylbutanal             |                              |         | 2.49  | 2.44    | 29, 41, 27, 57, 58, 39... 86    |
| 2-Butenal                   | 4.05                         | 4.00    |       |         | 41, 70, 39, 69, 29, 44          |
| 2-Methyl-2-butenal          | 4.67                         | 4.58    |       |         | 55, 29, 84, 27, 39              |
| 2-Heptenal                  |                              |         | 5.63  | 5.52    | 41, 83, 55, 57, 56, 39... 112   |
| 2-Octenal                   | 8.00                         | 7.95    |       |         | 41, 29, 27, 55, 39, 70... 126   |
| 2-Phenyl-2-butenal          |                              |         | 8.94  | 9.10    | 146, 117, 115, 91, 116, 39      |
| 5-Methyl-2-phenyl-2-hexenal |                              |         | 11.02 | 11.12   | 188, 91, 117, 115, 43, 39       |
| Benzaldehyde                |                              |         | 5.79  | 5.80    | 77, 105, 106, 51, 50            |
| Phenylacetaldehyde          |                              |         | 6.65  | 6.81    | 91, 92, 39, 65, 120             |
| Ketones                     |                              |         |       |         |                                 |
| Acetone                     | 1.00                         | 1.00    |       |         | 43, 58, 27, 26, 42, 29          |
| 2-Butanone                  |                              |         | 1.74  | 1.65    | 43, 29, 27, 72, 26, 57          |
| 2-Pentanone                 |                              |         | 2.60  | 2.70    | 42, 27, 86, 41, 39, 58          |
| 2-Heptanone                 |                              |         | 4.89  | 4.89    | 43, 58, 27, 71, 29, 41... 114   |
| 2-Octanone                  | 6.56                         | 6.61    |       |         | 43, 58, 41, 71, 59, 27... 128   |
| 2-Decanone                  |                              |         | 8.00  | 8.00    | 58, 43, 59, 71, 41... 156       |
| 2-Tridecanone               | 11.70                        | 11.58   |       |         | 58, 43, 59, 71, 41, 57... 198   |
| 2-Heptadecanone             |                              |         |       |         | 58, 43, 71, 85, 69, 83... 254   |
| 3-Pentanone                 |                              |         | 2.75  | 2.81    | 57, 29, 27, 86, 26, 42          |
| 3-Octanone                  | 6.24                         | 6.26    |       |         | 43, 57, 29, 71, 72, 99... 128   |
| 4-Methyl-2-pentanone        |                              |         | 3.47  | 3.28    | 43, 58, 57, 41, 29, 85... 100   |
| 3-Buten-2-one               | 2.94                         | 2.80    |       |         | 55, 43, 27, 70, 26, 42          |
| 3-Methyl-3-buten-2-one      |                              |         | 2.55  | 2.50    | 43, 41, 84, 69, 39, 42          |
| 3-Penten-2-one              | 4.90                         | 4.85    |       |         | 69, 41, 84, 43, 39, 40          |
| 2,3-Butanedione             |                              |         | 1.69  | 1.82    | 43, 86, 42, 29, 41              |
| 2,3-Pentanedione            |                              |         | 2.84  | 2.90    | 43, 29, 57, 27, 15, 100         |
| Cyclopentanone              | 5.55                         | 5.42    |       |         | 55, 28, 41, 84, 42, 56          |
| 2-Cyclopenten-1-one         | 7.30                         | 7.11    |       |         | 39, 82, 27, 54, 53, 26          |
| 1-Phenyl-2-propanone        |                              |         | 7.49  | 7.58    | 43, 91, 92, 134, 65, 39         |
| 1-Hydroxy-2-propanone       |                              |         |       | 3.37    | 43, 31, 74, 29, 42, 45          |

Table I (continued)

| Compound                                | Retention index <sup>a</sup> |         |                   |         | Mass spectral data <sup>b</sup> |
|-----------------------------------------|------------------------------|---------|-------------------|---------|---------------------------------|
|                                         | CBW                          |         | SF-96             |         |                                 |
|                                         | Known                        | Unknown | Known             | Unknown |                                 |
| 3-Hydroxy-2-butanone                    |                              |         | 3.08              | 3.50    | 45, 43, 27, 29, 18, 28...88     |
| 3-Hydroxy-2-pentanone                   |                              |         | 4.25              | 4.38    | 59, 31, 43, 45, 41, 58...102    |
| 2-Hydroxy-3-pentanone                   |                              |         | 4.27              | 4.45    | 45, 59, 31, 43, 29, 57...102    |
| Esters                                  |                              |         |                   |         |                                 |
| Methyl formate                          |                              |         | <1.00             | <1.00   | 31, 29, 32, 60, 15              |
| Ethyl formate                           | 1.00                         | 1.00    |                   |         | 31, 29, 27, 45, 74              |
| Ethyl acetate                           |                              |         | 2.00              | 2.00    | 43, 29, 27, 28, 45...88         |
| Pentyl acetate                          | 5.40                         | 5.35    |                   |         | 43, 70, 55, 15, 41, 27...130    |
| Acetol acetate                          |                              |         | 4.89              | 5.10    | 43, 29, 42, 116, 86, 27         |
| Ethyl propionate                        |                              |         | 3.00              | 3.00    | 29, 57, 27, 28, 26, 45...102    |
| Ethyl isovalerate                       | 4.28                         | 4.30    |                   |         | 29, 57, 27, 88, 41, 85...130    |
| Ethyl dodecanoate                       | 12.00                        | 12.00   |                   |         | 88, 29, 101, 41, 43, 27...228   |
| Ethyl tetradecanoate                    |                              |         | 14.00             | 14.00   | 88, 43, 41, 29, 101, 55...256   |
| Methyl hexadecanoate                    |                              |         |                   |         | 74, 87, 75, 69, 143, 83...270   |
| Acids                                   |                              |         |                   |         |                                 |
| Acetic acid                             |                              |         |                   |         | 43, 45, 60, 28, 15, 42          |
| Octadecanoic acid                       |                              |         |                   |         | 43, 57, 41, 55, 71, 73...284    |
| Lactones                                |                              |         |                   |         |                                 |
| $\delta$ -Nonalactone                   |                              |         | 10.32             |         | 99, 42, 27, 41, 29, 71...156    |
| $\delta$ -Decalactone                   | 15.28                        |         |                   |         | 99, 27, 42, 41, 71, 43...170    |
| Furans                                  |                              |         |                   |         |                                 |
| 2-Methylfuran                           | 2.00                         | 1.95    |                   |         | 82, 53, 81, 39, 27, 29          |
| 2-Ethylfuran                            |                              |         | 3.00              | 3.00    | 81, 39, 96, 41, 51, 65          |
| 2-Butylfuran                            |                              |         | 5.02              | 5.16    | 81, 82, 124, 53, 39, 41         |
| 2-Pentylfuran                           |                              |         | 6.00              | 5.85    | 81, 82, 138, 53, 41, 39         |
| 2-Furaldehyde                           |                              |         | 4.62              | 4.91    | 39, 96, 95, 29, 38              |
| 5-Methyl-2-furaldehyde                  |                              |         | 5.61              | 6.00    | 110, 109, 53, 27, 29, 51        |
| Furfuryl alcohol                        |                              |         | 5.33              | 5.90    | 39, 98, 41, 29, 81, 53          |
| 2-Acetylfuran                           |                              |         | 5.20              | 5.53    | 95, 110, 43, 39, 96, 68         |
| Propionylfuran                          |                              |         | 6.39              | 6.40    | 95, 124, 39, 96, 67, 55         |
| Furfuryl methyl ketone                  |                              |         | 5.67              | 5.80    | 43, 81, 53, 82, 124             |
| 2-Methyltetrahydrofuran-3-one           |                              |         | 4.32              | 4.18    | 43, 72, 100, 29, 45, 27         |
| Furfuryl formate                        |                              |         | 5.32              | 5.46    | 81, 53, 126, 52, 80, 27         |
| Furfuryl acetate                        |                              |         | 5.88              | 6.11    | 81, 43, 98, 52, 53, 140         |
| Furfuryl propionate                     |                              |         | 7.06              | 7.03    | 81, 27, 29, 52, 53, 98...154    |
| Furfuryl butyrate                       | 10.18                        | 9.81    |                   |         | 81, 27, 98, 43, 53, 52...168    |
| Furfuryl pentanoate                     |                              |         | 8.85              | 8.26    | 81, 27, 29, 98, 53, 41...182    |
| Furfuryl hexanoate                      | 12.00                        | 12.02   |                   |         | 81, 27, 98, 53, 43, 29...196    |
| Ethyl furoate                           | 9.79                         | 9.81    |                   |         | 95, 29, 39, 112...140           |
| Ethyl furfuryl ether                    |                              |         | 5.09              | 5.08    | 81, 82, 53, 27, 29, 126         |
| Furfuryl ether                          |                              |         | 9.38              | 9.42    | 81, 82, 53, 27, 97, 39...178    |
| 2,2'-Methylenedifuran                   |                              |         | 6.87              | 7.16    | 148, 91, 120, 147, 39, 65       |
| 5-Methyl-2,2'-methylenedifuran          |                              |         | 7.84              | 7.95    | 162, 91, 43, 147, 161, 119      |
| 2(or 3)-Phenylfuran                     |                              |         | 8.40 <sup>c</sup> | 8.64    | 144, 115, 145, 63, 116, 89      |
| Thiazoles                               |                              |         |                   |         |                                 |
| Thiazole                                | 6.18                         | 6.25    |                   |         | 58, 85, 57, 45, 26, 32          |
| 4-Methylthiazole                        | 6.50                         |         |                   |         | 99, 71, 72, 45, 39, 69          |
| 2-Acetylthiazole                        |                              |         | 6.30              | 6.62    | 43, 127, 99, 58, 57, 112        |
| Thiophenes                              |                              |         |                   |         |                                 |
| 2-Methylthiophene                       | 4.62                         | 4.58    |                   |         | 97, 98, 45, 39, 99, 27          |
| 2-Thiophenecarboxaldehyde               |                              |         | 6.19              | 6.54    | 111, 112, 39, 29, 45, 83        |
| 3-Thiophenecarboxaldehyde               |                              |         | 6.27              | 6.62    | 111, 112, 39, 45, 83, 57        |
| 3-Methyl-2-thiophene-carboxaldehyde     |                              |         | 7.48              | 7.21    | 125, 126, 97, 45, 53, 127       |
| 5-Methyl-2-thiophene-carboxaldehyde     |                              |         | 7.81              | 7.75    | 125, 126, 53, 97, 45, 127       |
| 2,5-Dimethyl-3-thiophene-carboxaldehyde |                              |         | 8.23              | 8.27    | 140, 139, 111, 59, 29, 45       |
| 2-Acetylthiophene                       |                              |         | 6.90              | 7.44    | 111, 126, 39, 43, 45, 83        |
| 3-Acetylthiophene                       |                              |         |                   |         | 111, 43, 126, 39                |
| 5-Methyl-2-acetylthiophene              |                              |         | 8.36              | 8.48    | 125, 140, 53, 45, 43, 97        |
| Thiopheneacrolein                       | 14.92                        | 15.52   |                   |         | 138, 45, 110, 29, 109, 39       |
| Sulfur Compounds                        |                              |         |                   |         |                                 |
| Methanethiol                            |                              |         | <1.00             | <1.00   | 47, 48, 45, 46, 15, 44          |
| Methyl sulfide                          | <1.00                        | <1.00   |                   |         | 62, 47, 45, 46, 61              |
| Methyl disulfide                        |                              |         | 3.30              | 3.47    | 94, 45, 79, 46, 47, 43          |
| Methyl sulfone                          | 12.12                        | 12.10   |                   |         | 79, 94, 45, 29, 48, 63          |
| Methyl sulfoxide                        | 9.29                         | 9.41    |                   |         | 63, 78, 45, 29, 46, 61          |
| Methyl thioacetate                      |                              |         | 3.00              | 3.00    | 43, 90, 45, 47, 48, 42          |
| Methyl thiopropionate                   | 4.93                         | 4.92    |                   |         | 29, 57, 104, 27, 45, 47         |
| Ethyl thiopropionate                    | 5.31                         | 5.30    |                   |         | 57, 29, 118, 27, 61, 58         |
| Benzyl methyl sulfide                   | 10.29                        | 10.26   |                   |         | 91, 138, 45, 65, 39, 92         |
| Furfuryl methyl disulfide               | 8.58                         | 8.38    |                   |         | 81, 128, 53, 27, 45, 51         |
| Furfuryl methyl disulfide               |                              |         | 8.39              | 8.46    | 81, 53, 27, 45, 51, 160         |

Table I (continued)

| Compound                                             | Retention index <sup>a</sup> |         |                   |         | Mass spectral data <sup>b</sup> |
|------------------------------------------------------|------------------------------|---------|-------------------|---------|---------------------------------|
|                                                      | CBW                          |         | SF-96             |         |                                 |
|                                                      | Known                        | Unknown | Known             | Unknown |                                 |
| Pyrroles                                             |                              |         |                   |         |                                 |
| Pyrrole-2-carboxaldehyde                             |                              |         | 8.44              | 8.27    | 95, 94, 39, 66, 38, 28          |
| 5-Methylpyrrole-2-carboxaldehyde                     |                              |         | 8.32              | 8.53    | 109, 108, 80, 53, 27, 29        |
| 1-Acetylpyrrole                                      | 9.06                         | 8.95    |                   |         | 67, 109, 39, 43, 41, 40         |
| 2-Acetylpyrrole                                      |                              |         | 8.39              | 8.18    | 94, 109, 39, 66, 38, 43         |
| 2-Propionylpyrrole                                   | 13.56                        | 13.70   |                   |         | 94, 123, 66, 39, 95, 28         |
| Pyrazines                                            |                              |         |                   |         |                                 |
| Pyrazine                                             |                              |         | 3.30              | 3.50    | 26, 80, 53, 28, 52, 51          |
| Methylpyrazine                                       |                              |         | 4.38              | 4.38    | 94, 67, 26, 39, 40, 55          |
| Ethylpyrazine                                        |                              |         | 5.21              | 5.30    | 107, 108, 26, 27, 39, 79        |
| 2,3-Dimethylpyrazine                                 |                              |         | 5.45              | 5.30    | 67, 108, 42, 41, 39, 52         |
| 2,5-Diethylpyrazine                                  |                              |         | 5.17              | 5.21    | 42, 108, 39, 38, 27, 81         |
| 2,6-Dimethylpyrazine                                 |                              |         | 5.12              | 5.28    | 108, 42, 39, 38, 67, 41         |
| 2-Ethyl-5-methylpyrazine                             |                              |         | 5.96              | 6.00    | 121, 122, 39, 27, 56            |
| 2-Ethyl-6-methylpyrazine                             |                              |         | 5.91              | 5.95    | 121, 122, 39, 42, 40, 27        |
| 2,5-Diethylpyrazine                                  |                              |         | 6.93              | 7.00    | 136, 121, 135, 39, 27           |
| Trimethylpyrazine                                    |                              |         | 6.16              | 6.09    | 42, 122, 39, 27, 81, 54         |
| 2-Ethyl-3,5-dimethylpyrazine                         |                              |         | 6.83              | 6.93    | 135, 136, 39, 42, 56, 54        |
| 2-Ethyl-3,6-dimethylpyrazine                         |                              |         | 6.90              | 6.81    | 135, 136, 42, 56, 39, 108       |
| 5-Ethyl-2,3-dimethylpyrazine                         |                              |         |                   |         | 135, 136, 42, 39, 56, 108       |
| 2,3-Diethyl-5-methylpyrazine                         |                              |         | 7.63              | 7.59    | 150, 135, 56, 39, 149, 41       |
| 2,5-Diethyl-3-methylpyrazine                         | 8.70                         | 8.52    |                   |         | 149, 150, 39, 135               |
| 3,5-Diethyl-2-methylpyrazine                         | 8.60                         | 8.47    |                   |         | 149, 150, 39, 122, 53, 135      |
| Triethylpyrazine                                     | 9.06                         | 9.08    |                   |         | 164, 149, 163, 39, 56, 136      |
| Tetramethylpyrazine                                  | 8.41                         | 8.30    |                   |         | 54, 136, 42, 27, 52, 39         |
| 2-Methyl-3-vinylpyrazine                             | 6.28                         | 6.42    |                   |         | 119, 120, 67, 26, 39, 52        |
| 2-Methyl-5-vinylpyrazine                             |                              |         | 6.26              | 6.40    | 120, 52, 54, 39, 51, 27         |
| Acetylpyrazine                                       |                              |         | 6.30              | 6.66    | 43, 122, 80, 79, 52, 53         |
| 2-Methyl-5(or 6)-acetylpyrazine                      |                              |         | 7.27 <sup>d</sup> | 7.50    | 136, 43, 94, 93, 39, 67         |
| 2-Ethyl-3-acetylpyrazine                             |                              |         | 7.81              | 7.93    | 43, 107, 150, 52, 27, 79        |
| 2-Ethyl-5(or 6)-acetylpyrazine                       |                              |         | 7.70 <sup>e</sup> | 7.90    | 150, 43, 107, 108, 53, 27       |
| 6,7-Dihydro-5H-cyclopentapyrazine                    |                              |         | 7.18              | 7.43    | 120, 119, 39, 65, 52, 27        |
| 2-Methyl-6,7-dihydro-5H-cyclopentapyrazine           | 10.73                        | 10.60   |                   |         | 134, 133, 39, 66, 40, 107       |
| 5-Methyl-6,7-dihydro-5H-cyclopentapyrazine           |                              |         | 7.47              | 7.56    | 119, 134, 133, 39, 27, 52       |
| 2(or 3)-5-Dimethyl-6,7-dihydro-5H-cyclopentapyrazine | 10.26                        | 10.26   |                   |         | 148, 66, 147, 39, 27, 107       |
| Quinoxaline                                          | 12.71                        | 12.55   |                   |         | 130, 76, 103, 50, 75, 51        |
| 2-Methylquinoxaline                                  |                              |         | 9.24              | 9.39    | 144, 117, 76, 50, 77            |
| 6-Methylquinoxaline                                  |                              | 13.08   |                   |         | 144, 90, 89, 117, 143, 145      |
| 5,6,7,8-Tetrahydroquinoxaline                        |                              |         | 8.05              | 8.22    | 134, 133, 52, 106, 39, 27       |
| 2-Methyl-5,6,7,8-tetrahydroquinoxaline               |                              |         | 9.00              | 9.05    | 148, 147, 39, 52, 120, 79       |
| (2-Furyl)pyrazine                                    | 12.80                        | 13.35   |                   |         | 146, 93, 63, 64, 92, 39         |
| 2-(2-Furyl)-5(or 6)-methylpyrazine                   |                              |         | 9.87              | 9.64    | 160, 92, 63, 39, 93, 64         |
| Miscellaneous                                        |                              |         |                   |         |                                 |
| Ethyl vinyl ether                                    |                              |         | <1.00             | <1.00   | 44, 43, 72, 29, 27              |
| Diethyl ether                                        |                              |         | 1.10              | 0.97    | 31, 59, 29, 45, 27, 74          |
| Acetaldehyde, diethyl acetal                         |                              |         | 3.20              | 3.18    | 45, 73, 29, 27, 43...118        |
| Phenol                                               | 13.46                        | 13.40   |                   |         | 94, 39, 66, 65, 40              |
| Trimethyloxazole                                     | 5.50                         | 5.68    |                   |         | 43, 42, 111, 55, 68, 27         |
| Vanillonitrile                                       |                              |         |                   |         | 134, 149, 106, 29, 51, 78       |

<sup>a</sup> Compounds with an  $I_E$  value lower than that of ethyl formate (1.00) are given as <1.00. <sup>b</sup> These data were taken from the mass spectra of known compounds. However, in every case the data for the unknown compound matched that for the known. The mass spectrometer used had a source temperature of 150°, an ionizing voltage of 70 eV, and an ionizing current of approximately 3.5 A. Numbers in italics are the molecular weights of the compounds. <sup>c</sup> The 2 and 3 isomers have very similar retention indices and mass spectra. The authentic compounds were synthesized by E. Klaiber and A. O. Pittet of IFF. <sup>d</sup>  $I_E$  for 2-methyl-5-acetylpyrazine.

and trap 6 was analyzed on a 3 ft  $\times$   $\frac{1}{8}$  in. stainless steel column packed with 2% OV-17 on 80-100 mesh acid-washed, DMCS-treated Chromosorb W. All capillary analyses were carried out at a helium flow rate of 10 ml/min. The oven temperature was programmed from 70 to 190° at 1°/min for the total extract and trap 1, and at 2°/min for traps 2-5.

Identifications were based on the comparison of known and unknown mass spectra and confirmed wherever possible by determining the retention indices, or  $I_E$  values, relative to a series of ethyl esters of normal alkanic acids (van den Dool and Kratz, 1963). Additional confirmation for certain compounds was provided by the gc analysis of the total extract and each trap for the presence of sulfur and/or nitrogen containing compounds. The instrument used was a Tracor MT 220 gas chromatograph modified in our laboratory as described by Walradt (1973).

For identification purposes, some compounds for which mass spectral data were not available were synthesized: furfuryl methyl ketone (Hass *et al.*, 1950); ethyl furfuryl ether (Kirner, 1928); propionylfuran (Heid and Levine, 1948); furfuryl methyl disulfide (Milligan and Swan, 1963); 2-thiophene acrolein (Hori, 1958); 2,4,5-trimethyl-oxazole (Reppe and Magin, 1963); ethyl thiopropionate, by a modification of the procedure of McFadden *et al.* (1965); and furfuryl propionate, by a modification of a method described by Blatt (1941). The bicyclic pyrazines were synthesized by Pittet *et al.* (1974).

## RESULTS AND DISCUSSION

Table I lists the compounds that were identified. In some cases, the small amount of sample available precluded the determination of  $I_E$  values. In other cases the "known" compound was not readily available. In the case of some very polar compounds, for example furfuryl alcohol, the known and unknown  $I_E$  values may differ considerably due to the fact that the  $I_E$  value of the unknown is determined on a complex mixture "spiked" with ethyl esters where its retention time may be affected by the other components. The  $I_E$  value for the known compounds, however, is determined on a mixture of the pure compound with a series of ethyl esters. This effect is particularly evident with free acids. Where there is a discrepancy in the  $I_E$  value, the compound was only reported if its mass spectral fragmentation pattern was very characteristic. In those cases where the position of the double bond was not specified, this information could not be unambiguously determined from the data.

Table II shows some quantitative data for the major classes of compounds found in cooked pork liver. The figures were determined by computer area normalization and incorporate data for both the positively identified com-

**Table II. Relative Per Cent of Major Classes of Compounds in Volatiles from Cooked Pork Liver**

| Class            | Quantity, % |
|------------------|-------------|
| Ketones          | 0.53        |
| Esters           | 0.56        |
| Hydrocarbons     | 1.33        |
| Miscellaneous    | 1.38        |
| Sulfur compounds | 2.16        |
| Alcohols         | 4.21        |
| Aldehydes        | 14.07       |
| Furans           | 28.80       |
| Pyrazines        | 40.96       |
| Total            | 94.00       |

pounds reported in this paper as well as some tentatively identified compounds. It should be pointed out that pyrazines were by far the largest class of compounds found, accounting for almost 41% of the isolated components amenable to gc analysis.

## SUMMARY

Pork liver, cooked superatmospherically at 325°F, was analyzed by gc and coupled gc-mass spectrometry. A total of 179 components have been identified. Although none of these chemicals can be considered to have a typical liver character, many of them are believed to contribute to the overall flavor quality.

## ACKNOWLEDGMENT

The authors thank Anne Sanderson for assistance in obtaining and interpreting mass spectra.

## LITERATURE CITED

- Blatt, A. H., Ed., "Organic Syntheses," Collect. Vol. I, 2nd ed., Wiley, New York, N. Y., 1941, p 285.  
 Hass, H. B., Susie, A. G., Heider, R. L., *J. Org. Chem.* **15**, 8 (1950).  
 Heid, J. V., Levine, R., *J. Org. Chem.* **13**, 409 (1948).  
 Hori, I., *J. Sci. Res. Inst. Tokyo* **52** (1479), 78 (1958).  
 Kirner, W. R., *J. Amer. Chem. Soc.* **50**, 1955 (1928).  
 McFadden, W. H., Seiffert, R. M., Wasserman, J., *Anal. Chem.* **37**, 560 (1965).  
 Milligan, B., Swan, J. M., *J. Chem. Soc. A*, 6008 (1963).  
 Pittet, A. O., Muralidhara, R., Walradt, J. P., Kinlin, T. E., *J. Agr. Food Chem.*, **22**, 273 (1974).  
 Reppe, W., Magin, A., French Patent 1,340,996 (Sept 16, 1963).  
 van den Dool, H., Kratz, P. D., *J. Chromatogr.* **11**, 463 (1963).  
 Walradt, J. P., in "Techniques of Combined Gas Chromatography/Mass Spectrometry: Applications in Organic Analysis," McFadden, W. H., Ed., Wiley, New York, N. Y., 1973.  
 Watson, J. T., Biemann, K., *Anal. Chem.* **37**, 844 (1965).  
 Williams, A. A., *Chem. Ind. (London)*, 1510 (1969).

Received for review February 13, 1974. Accepted May 6, 1974. Presented at the 166th National Meeting of the American Chemical Society, Chicago, Ill., Aug 26-31, 1973.