

РП: S0960-894X(96)00452-0

THE SYNTHESIS OF 5-AZIDO-UTP AND UDP-PYRIDOXAL AND THEIR EFFECT ON THE CARBAMOYL-PHOSPHATE SYNTHETASE (CPS II) IN THE MAMMALIAN MULTIENZYME POLYPEPTIDE CAD

Li-Min Zhu * a and Elizabeth A. Carrey *

Department of Biochemistry, University of Dundee, Dundee DD1 4HN, Scotland, U.K.

Abstract: 5-azidouridine-5'-triphosphate and UDP-pyridoxal have been synthesized using modified synthetic methods. The two synthetic UTP analogues, with commercially available analogues and chemicals have been tested for their effect on CPSase II in the mammalian multienzyme polypeptide CAD. Structure activity relationships are discussed.


Copyright © 1996 Elsevier Science Ltd

Pyrimidine biosynthesis in mammalian cells is initiated by the three enzyme activities found in the multifunctional polypeptide CAD [reviewed in 1]. The three enzymes are the glutamine-dependent carbamoyl-phosphate synthetase (CPSaseII), aspartate transcarbamoylase (ATCase) and dihydroorotase (DHOase). CPSaseII is allosterically controlled by UTP (negative effector) and phosphoribosyl pyrophosphate PRibPP (positive effector) [2,6].

Reactive nucleotide analogues have provided valuable information on the structure and mechanism of a number of enzymes. The mechanism by which the interaction of metabolites with the enzyme leads to modulation of its activity are clearly of interest. Furthermore, the biosynthesis of pyrimidines is essential for most growing cells, thus the design and synthesis of analogues for the metabolites are a potential target for the drug discovery for chemotherapy against cancer. We have used two synthetic analogues, other nucleotides and chemicals to compare the structures of effectors and the enzyme activity of the glutamine-dependent carbamoyl-phosphate synthetase (CPSaseII) domain of CAD, to examine the effect of UTP analogues on CPSase II in the multienzyme polypeptide CAD and demonstrate the importance of the structures of effector analogues to the CPSase II activities. The studies will increase our understanding of the binding sites of PRibPP and UTP in the mammalian multienzyme polypeptide CAD, which regulate the biosynthesis of pyrimidines [3].


The synthetic routes to UDP-pyridoxal are outlined in Scheme 1. UDP-pyridoxal **3** was synthesized by a modification of previously published methods [4a, 4b] and from the coupling of pyridoxal 5-phosphate **1** and UMP-morpholidate **2** in anhydrous pyridine. The extracted aqueous layer was applied to an anion exchange column of Dowex 1-X8 (Cl⁻ form) to remove UMP and pyridoxal 5-phosphate, the yield was 65%. 5-azidouridine 5'-triphosphate

^a Present address: School of Molecular and Life Sciences, University of Abertay Dundee, Dundee DD1 1HG. Fax: +44 1382 308663, E-mail: mltlz@tay.ac.uk.



Reagents and conditions: a) Pyridine (anhydrous), 3-5 days, rt; b) Extracted with ethylether; c) Dowex 1-X8 (Cl⁻ form), washed with H₂O, 3mM-10mM LiCl, eluted with 10mM HCl; d) Adjusted pH to 6.5 with 1N LiOH, kept 0°C, 1h, centrifuged, dried.

Scheme 2




Reagents and conditions: a) Dowex 50W resin column (H^+ form); b) NOBF₄/DMF (anhydrous); c) 1N HCl, 15 min, neutralized by NH₄OH; d) Sephadex G-10 column (H₂O); e) Benzyl diethylaminoethyl cellulose column (gradient of NH₄HCO₃); f) Zn/HCl, 0.5-2h, neutralized by NH₄OH; g) Diethylaminoethyl cellulose column (gradient of NH₄HCO₃); h) Sephadex G-10 column (H₂O); i) NaN₃, 10 min, 0°C, 30 min, rt; j) Sephadex G-10 column (H₂O); k) Benzyl diethylaminoethyl cellulose column (gradient of NH₄HCO₃); l) Dowex 50W column (H₂O). Overall yield:37.5% [5].

(N₃-UTP) 7 was chemically synthesized by modification of a previously published synthetic method for 5-azido-2'-deoxyuridine 5'-mono-phosphate [4c] and the other chemical synthetic methods for nucleotides [4d-e]. The synthetic routes to N₃-UTP 7 are outlined in Scheme 2. UTP 4 (sodium salt) was firstly converted to the free acid form by applying the nucleotide to a Dowex 50W resin (H⁺ form) column. The dried UTP was nitrated by nitrosonium tetrafluoroborate. The neutralized sample was desalted using Sephadex G-10 column and then an anion and cation mixed bed exchange benzyl diethylaminoethyl cellulose column to remove mixed ions. 5-NO₂UTP 5 was reduced to 5-NH₂UTP 6 by adding the mixture of clean granulated zinc and HCl. Then the neutralized 5-NH₂UTP 6 was applied to an anion exchange diethylaminoethyl cellulose column and desalted by Sephadex G-10 gel filtration chromatography. 5-NH₂UTP 6 was finally converted to the product N₃-UTP 7 by the addition of NaN₃ which was produced from the mixture of NaNO₂ and H₂O. The neutralized sample was again desalted by gel filtration chromatography on a Sephadex G-10 column, a mixed bed benzyl diethylaminoethyl cellulose column and a strongly acidic cation Dowex 50W column to remove remove residual NH₄⁺. The overall yield was 37.5% from UTP.

The two synthetic UTP analogues UDP-pyridoxal 3 and N₃-UTP 7, with commercially available analogues and chemicals: uridine 5'-diphosphate (UDP) 8, uridine (Urd) 9, pyridoxal 5-phosphate 12, pyridoxal hydrochloride 13, D-Ribose-5-phosphate 11 and 1'-phospho-ribosyl-5'-pyrophosphate (PRibPP) 10 have been tested for their effect on CPSase II in CAD. Table I shows the results that 5-azidouridine-5'-triphosphate 7, uridine 5'diphosphate 8 and pyridoxal hydrochloride 13, like UTP 1, inhibit the activity of the glutamine-dependent carbamyl phosphate synthetase (CPSase II) in the mammalian multienzyme polypeptide CAD. However, UDP-pyridoxal 3, pyridoxal 5-phosphate 12, like PRibPP 10, activate CPSase II, while uridine 9 and D-ribose 5-phosphate 11 have no effect. Tamura et al. have demonstrated that ADP-pyridoxal is a potent inhibitor of yeast alcohol dehydrogenase and rabbit muscle pyruvate kinase [4b], which promted our interest to use UDP-pyridoxal as an UTP analogue to inhibit the activity of CPS II. However, the results show that UDP-pyridoxal, not as expected as an inhibitor, activates the enzyme. Comparing the structures of UTP 1 and UDP-pyridoxal 3, the major difference between two compounds is that pyridoxal group in the UDP-pyridoxal replaces the γ -position phosphate in UTP, to cause the totally opposite effect on the CPSase II activities. Obviously, the structure of this part of the effectors is really important to the activity of the enzyme. Table I demonstrates that UDP 8 is also a good inhibitor of the activity of the enzyme even though lacking one phosphate group compared to UTP 1. However, very interestingly, Uridine 9 has completely no effect, which suggests that phosphate groups in UTP are the most important groups to inhibit the activity on CPSase II in the mammalian multienzyme polypeptide CAD. On the other hand, where the 5-H in UTP 1 has been replaced by an azido group, the UTP analogue N_3 -UTP 7, inhibits the activity of CPSase II, but is much less efficiency than UTP. Table I shows that 52% activity of the CPSase II has been inhibited by the same concentration of N3-UTP 7 and 82% by UTP 1. The results suggest that the 5-H position or the base in UTP is one of the

No. of Entry	Nucleotides and chemicals		Activity	Concentration
1	ATP (Substrate only)		% [100]	mM [4]
2	+ 1		18	1.2
3	+ 7		48	1.2
4	+ 8	HO-P-O-CH, OH OH OH	54	1.2
5	+ 13	HO HO H,C N HCI	73	1.2
6	+ 9	HOCH ₂ O OH OH	100	1.6
7	+ 11	но но он он	100	1.6

 Table I. Effect of nucleotides and chemicals on the activity of CPSase II in the multienzyme polypeptide CAD



The activities were determined using the method given in [7]. pH for all used analogues and chemicals were adjusted to 7.5 before use.

most important positions to cause the inhibition of the enzyme activity. Thus, a single change in structure at the both 5-H position and phosphate position of UTP has a dramatic effect, changing the enzyme activity.

We believe that the reason that UDP-pyridoxal 3 activates the enzyme is the presence of pyridoxal 5-phosphate 12. Table I shows that pyridoxal 5-phosphate has increased the activity of the CPSase II by 70% and UDP-pyridoxal by only 20%. The results suggest that the partial structure of UDP-pyridoxal 3, the pyridoxal 5-phosphate 12, is the cause for the activation of the enzyme activity. Pyridoxal hydrochloride 13 has been used to investigate if the pyridoxal group activates the enzyme activity, and the result demonstrates that pyridoxal itself, without one more crucial phosphate group, inhibits the activity of the enzyme (27% inhibition). The evidences suggest, once again, that the phosphate group of UTP is the most important group to act as an efficient allosteric inhibitor. The phosphate group has a crucial role not just for the allosteric inhibitor UTP, but also for the allosteric activator PRibPP, to regulate the biosynthesis of pyrimidine. The experimental results demonstrate that the analogue D-Ribose-5-phosphate 11, unlike PRibPP 10, has no effect on the CPSase II. Comparing the structures of both PRibPP 10 and D-Ribose-5-phosphate 11, it is not difficult to draw a conclusion that D-Ribose-5-phosphate does not inhibit or activate the enzyme because D-Ribose-5-phosphate has two missing phosphate groups, of which one is 1'-phosphate and another is 5'-phosphate in PRibPP.

The two synthetic analogues UDP-pyridoxal 3 and N_3 -UTP 7, with commercially available analogues and chemicals used on the activity of CPSase demonstrate how important the structure of effectors is to the activity of the enzyme. The results will certainly help the design of more efficient analogues for the metabolites which regulate the biosynthesis of pyrimidine and the drug discovery for chemotherapy against cancer.

Acknowledgements: this work was supported by the Medical Research Council. We are grateful to: Dr Robert A. Field, School of Chemistry, University of St Andrews; Professor J.S. Brimacombe, Department of Chemistry, University of Dundee.

References and Notes

- (a) Jones, M.E. Annu. Rev. Biochem. 1980, 49, 253; (b) Carrey, E.A. Biochem. Soc. Trans. 1993, 21, 191; (c) Davidson, J.N., Chen, K.C., Jamison, R.S., L.A., Musmanno, L.A. & Kern, C.B. BioEssays, 1993, 3, 157; (d) Evans, D.R., Bein, K., Guy, H.I., Liu, X., Molina, J.A. and B.H. Zimmermann, B.H. Biochem. Soc. Trans. 1993, 21, 186.
- (a) Evans, D.R. Multidomain Proteins: Structure and Evolution (Hardie, D.G. & Coggins, J.R., eds) Elsevier Biomedical Press, Amsterdam, 1986; (b) Carrey, E.A. D.G. Hardie, D.G. Eur. J. Biochem., 1988, 171, 583; (c) Carrey, E.A., Campbell, D.G. & Hardie, D.G. EMBO J. 1985, 4, 3735; (d) Carrey, E.A. Biochem. J. 1986, 236, 327; (e) Hemmens, B. & Carrey, E.A. Eur. J. Biochem. 1994, 225, 845.
- (a) Devaney, M.A. & Powers-Lee, S.G. J. Biol. Chem. 1984, 259, 703; (b) Zhu, L.-M. & Carrey, E.A. Biochem. Soc. Trans. 1995, 23, 620.
- (a) Tagaya, M., Nakano, K. & Fukui, T. J. Biol. Chem. 1985, 260, 6670; (b) Tamura, J.K., Rakov, R.D. & Cross, R.L. J. Biol. Chem. 1986, 261, 4126; (c) Evans, R.K. & Haley B.E. Biochem. 1987, 26, 269; (d) Michelson, A.M. Biochim. Biophys. Acta. 1964, 91, 1; (e) Boulay, F., Dalbon, P. & Vignais, P.V. Biochem. 1985, 24, 7372.
- Abbreviations: CAD, the multienzyme polypeptide carrying the activities carbamoyl phosphate synthetase II aspartate transcarbamoylase and dihydroorotase; BD-Cellulose, benzyl diethylaminoethyl-Cellulose; DEAE-Cellulose, diethylaminoethyl-Cellulose; 5-N₃UTP, 5-Azidouridine-5'-triphosphate; PRibPP, 1'-phospho-ribosyl-5'-pyrophosphate.
- Enzymes: aspartate transcarbamoylase (EC 2.1.3.2., ATCase); glutamine-dependent carbamoyl-phosphate synthase (EC 6.3.4.16., CPSase II); dihydroorotase (EC 3.5.2.3. DHOase). the baby hamster kidney cell line 165-28 was a gift from Professor D.R. Evans, Wayne State University, Detroit.
- 7. Assay procedure: Carbamoyl-phosphate synthase was assayed by a modification of previously methods [2]. A total incubation volume of 250µl contained 0.1M Tris, pH 7.5 at 37°C, 0.1M KCl, 7.5% dimethyl sulphoxide (DMSO), 2.5% glycerol, 5 mM MgCl₂, 4 mM ATP, 1 mM dithiothreitol, and nitrogen donor either glutamine (3 mM) or NH₄Cl, pH 7.5 (40 mM). The reaction was initiated by adding 10µl of CAD solution (1-3µg of protein) and 15 µl of NaH¹⁴CO₃ (0.2M, 0.75Ci/mol). After 15 min, the assay was stopped by mixing with 0.10 ml of 2 M NH₄Cl and heating the tubes at 90°C for 20 min. The unchanged NaHCO₃ was driven off as CO₂ by adding 0.1 ml of 14% (v/v) HClO₄, and the assay mixtures were dried and resuspended for scintillation counting.

(Received in Belgium 9 July 1996; accepted 23 September 1996)