

Efficient and Selective Epoxidation of Alkenes with Sodium Periodate Using Supported Manganese Porphyrins under Ultrasonic Irradiation

Shahram Tangestaninejad* and Valiollah Mirkhani
Chemistry Department, Esfahan University, Esfahan 81744, Iran

(Received August 14, 1998; CL-980625)

Epoxides are readily obtained in high yields and good selectivities by ultrasonic irradiation of alkenes and sodium periodate in the presence of catalytic amounts of manganese porphyrins supported on polyvinylpyridine and IRA-900 ion-exchange resin.

Synthetic metalloporphyrins, like cytochrome P-450, act as useful catalysts for a variety of redox reactions.¹⁻³ The high efficiency of some of these catalytic systems makes them potentially useful for large scale oxidations. Immobilization of expensive metalloporphyrin catalysts on insoluble organic and inorganic supports appear to be a good way to render them practicable and improve their stability and selectivity.⁴⁻⁹ However, the activity of such heterogeneous catalysts is often less than that of soluble ones, due to either chemical or diffusional restrictions. The successful applications of ultrasonic irradiation in a number of heterogeneous reactions,^{10,11} prompted us to investigate the effect of ultrasonic irradiation on oxygenation reactions catalyzed by supported metalloporphyrin catalysts which are usually easy to handle but less reactive.

In this communication, we report a significant improvement in efficiency of supported metalloporphyrin catalysts through the use of ultrasonic irradiation. Higher catalytic activity was observed compared to the previously described homogeneous manganese porphyrin-periodate catalytic systems.^{12,13}

Using cyclooctene as the typical alkene substrate, the two catalysts i.e. sulfonated manganese(III) tetraphenylporphyrin supported on Amberlite IRA-900 ion exchange resin (MnTPPS-Ad)¹⁴ and sulfonated manganese(III) tetraphenylporphyrin supported on poly(4-vinylpyridine) (MnTPPS-PVP)¹⁵ were examined under various reaction conditions. The results are summarized in Table 1. We found that sonication has general applicability for supported metalloporphyrin catalysts with different oxygen donors. Higher activity and stability of the supported catalysts was observed with sodium periodate compared to other oxygen donors under ultrasonic irradiation.

The effect of different common organic solvents such as CH₃CN, (CH₃)₂CO, CH₃OH, CH₃CN/H₂O and (CH₃)₂CO/H₂O on epoxidation of cyclooctene with sodium periodate by supported manganese porphyrins was also studied. The same rate enhancement was obtained in all of these solvents under ultrasound irradiation. Among these solvents, CH₃CN/H₂O (1/1) mixture, which showed higher epoxide yield, was chosen as the reaction medium.

As shown in Table 2, sonication of a mixture of supported manganese porphyrins, alkene and NaIO₄ converts different alkenes efficiently to the corresponding epoxides.¹⁶ Typically, 1 mmole alkene, 2 mmole NaIO₄, 0.2 mmole imidazole, 0.0147 mmole supported manganese porphyrin in 20 mL CH₃CN/H₂O

Table 1. Catalytic epoxidation of cyclooctene by supported manganese porphyrins with different oxygen donors^a

Run	Catalyst	Oxygen donor	Time/min	Epoxide / % ^b	
				U ^c	S ^c
1	MnTPPS-PVP ^d	KHSO ₅	30	93	50
2	MnTPPS-Ad ^d	KHSO ₅	30	100	56
3	MnTPPS-PVP ^d	NaOCl	60	82	43
4	MnTPPS-Ad ^d	NaOCl	60	94	52
5	MnTPPS-PVP	NaIO ₄	60	100	35
6	MnTPPS-Ad	NaIO ₄	60	95	34

^aReaction conditions: cyclooctene (1 mmole), oxygen donor (2 mmole), imidazole (0.2 mmole), catalyst (mass equiv. to 0.0147 mmole Mn), CH₃CN / H₂O (10 mL/10 mL). ^bGLC yield based on starting cyclooctene. ^cU: ultrasonic irradiation (40 KHz), S: stirred magnetically (1200 rpm). ^dInstead of 10 mL H₂O, 10 mL of phosphate buffer at pH=7 was added.

(1/1) were sonicated in an ultrasound cleaner (40 KHz) at room temperature. Reactions were monitored by GLC and the identity of pure epoxides confirmed by ¹H NMR and MS spectral data.

Table 2. Epoxidation of alkenes with periodate using supported manganese porphyrins^a

Run	Alkene	Catalyst	Conversion/% ^b		Epoxide / % ^b	
			U ^c	S ^c	U ^c	S ^c
1	Cyclohexene	MnTPPS-Ad	92	40	88	30
2	Cyclohexene	MnTPPS-PVP	100	42	96	34
3	Styrene	MnTPPS-Ad	90	35	89	28
4	Styrene	MnTPPS-PVP	100	40	94	30
5	α -Methylstyrene	MnTPPS-Ad	100	45	92	38
6	α -Methylstyrene	MnTPPS-PVP	100	48	94	39
7	(+)-Camphene	MnTPPS-Ad	74	30	74	25
8	(+)-Camphene	MnTPPS-PVP	97	40	97	37
9	R-(+)-Limonene	MnTPPS-Ad	62	26	42 ^d	15 ^d
					20 ^e	11 ^e
10	R-(+)-Limonene	MnTPPS-PVP	84	34	57 ^d	24 ^d
					27 ^e	10 ^e
11	<i>trans</i> -Stilbene	MnTPPS-Ad	44	16	44 ^f	16 ^f
12	<i>trans</i> -Stilbene	MnTPPS-PVP	50	18	50 ^f	18 ^f
13	<i>cis</i> -Stilbene	MnTPPS-Ad	76	32	72 ^g	30 ^g
					4 ^f	2 ^f
14	<i>cis</i> -Stilbene	MnTPPS-PVP	80	33	76 ^g	30 ^g
					4 ^f	3 ^f

^aReaction conditions: alkene (1 mmole), NaIO₄ (2 mmole), imidazole (0.2 mmole), catalyst (mass equiv. to 0.0147 mmole Mn), CH₃CN/H₂O (10 mL/10 mL), time (1 h). ^bGLC yield based on starting alkene. ^cU: Ultrasonic irradiation (40 KHz), S: Stirred magnetically (1200 rpm). ^d1,2-Epoxide; ^e8,9-Epoxide; ^f*trans*-Epoxide; ^g*cis*-Epoxide.

To assess long-term stability for MnTPPS-PVP and MnTPPS-Ad catalysts under ultrasonic irradiation, recycling experiments were carried out with a single sample of each catalyst. Each of the catalysts were used repeatedly over styrene reactions, with the catalysts being removed by filtration at the end of each reaction, washed with water and acetonitrile and reused. The yields and leaching data obtained are listed in Table 3.

Table 3. Recycling data for supported manganese porphyrin catalysts in the epoxidation of styrene under ultrasonic irradiation^a

Run	Conversion / % ^b	Epoxide / % ^b	Mn leached % ^c
1	100(90) ^d	94(89)	2.2(2.5)
2	95(88)	89(87)	2.1(2.3)
3	91(85)	86(82)	1.2(1.3)
4	88(83)	82(80)	1.2(1.2)

^a Reaction Conditions: identical with those in Table 2, footnote a. ^b GLC yield based on starting styrene. ^c Reaction solutions assayed for leached Mn using atomic absorption spectroscopy and expressed as percentage of Mn originally present on resin. ^d The values in paranthesis are for MnTPPS-Ad / NaIO₄ and other values are for MnTPPS-PVP / NaIO₄ systems.

Alkene oxidations catalyzed by supported manganese porphyrins under ultrasound have two main advantages over those carried out under agitation with magnetic stirrer. The first is that oxidations are faster, typically complete within 1 h compared to at least 4 h required for agitation. Secondly, generally observed side reactions such as double bond cleavage and allylic oxidation are minimized. Further work concerning the extension of this reaction to various metalloporphyrin catalytic systems is currently under way.

This work was supported by special Grant(No. 751038) from University of Esfahan.

References and Notes

- 1 I. Tabushi, *Coord. Chem. Rev.*, **86**, 1 (1988).
- 2 B. Meunier, *Chem. Rev.*, **92**, 1411 (1992).
- 3 M. Sono, M. P. Roach, E. D. Coulter, and J. H. Dawson, *Chem. Rev.*, **96**, 2841 (1996).
- 4 D. R. Leonard and J. R. Lindsay Smith, *J. Chem. Soc., Perkin Trans. 2*, 25 (1991).
- 5 T. Mori, T. Santa, and M. Hirobe, *Tetrahedron Lett.*, **26**, 5555, (1985).
- 6 I. Mochida, K. Tsuji, K. Suetsugu, H. Fujitsu, and K. Takeshita, *J. Phys. Chem.*, **84**, 3159 (1980).
- 7 M. Nakamura, T. Tatsumi, and H. Tominage, *Bull. Chem. Soc. Jpn.*, **63**, 3334(1990).
- 8 S. Tangestaninejad and M. Moghadam, *Synth. Commun.*, **28**, 427 (1998).
- 9 S. Tangestaninejad and M. Moghadam, *J. Chem. Res. (S)*, **1998**, 242.
- 10 T. J. Mason, "Advances in Sonochemistry," JAI Press Ltd., London (1993), Vol. 3.
- 11 J. Ichihara and T. Hanafusa, *J. Chem. Soc., Chem. Commun.*, **1989**, 1843.
- 12 D. Mohajer and S. Tangestaninejad, *J. Chem. Soc., Chem. Commun.*, **1993**, 240.
- 13 D. Mohajer and S. Tangestaninejad, *Tetrahedron Lett.*, **35**, 945 (1994).
- 14 Y. Saito, M. Mifume, T. Kawaguchi, J. Odo, Y. Tanaka, M. Chikuma, and H. Tanaka, *Chem. Pharm. Bull.*, **34**, 2885 (1986).
- 15 S. Campestrini and B. Meunier, *Inorg. Chem.*, **31**, 1999 (1992).
- 16 In most epoxidation reactions under sonication or simple stirring the MnTPPS-PVP showed higher activity than MnTPPS-Ad catalyst. This suggests that the poly(vinylpyridine) resin support leads to a better activity for the catalyst compared to the simple cationic resin (Amberlite). The reason for this behavior is proximal effect due to coordination of a pyridine unit arising from the polymer.